

# Medical device

## Vascular Access System

AP123 V1.0

# Vascular Access System



## High Performance

- Designed for unparalleled support, balanced with flexibility based on its **unique LCT design**

## Procedural Flexibility

- New levels of procedural **efficiencies due to a unique ID/OD ratio**
- Large bore **ID (.087")** accommodates **.072"** aspiration catheters

## Maximized Control

- Confident **transmission accuracy** from tip to tail

# Elevate Performance with LCT

## 1 High-performance design

LCT technology allows for **infinite variations** on cut patterns of a single hypo tube. This **enables smoother transitions** and larger IDs **without increasing ODs** with a lower wall thickness. With **more metal and less plastic**, this design maximizes the catheter performance.

## 2 Procedural flexibility

The ability to produce incredibly thin walls allows for a larger inner lumen balanced with a reduced outer profile... the result, **a versatile choice that accommodates more devices than current technology on the market.**

## 3 Maximized control

Achieve **excellent pushability with great responsiveness** due to the customized cut pattern along the shaft without **sacrificing column strength or trackability.**

Tests performed by and data on file at Company .

# A **vascular** access system that provides **unparallel support**.



|                                                                                     |                                                                                    |                                           |                                          |                     |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------|------------------------------------------|---------------------|
|   | French Designation (FD)<br><b>6F</b>                                               | Outer Diameter (OD)<br><b>0.083in</b>     | Internal Diameter (ID)<br><b>0.074in</b> |                     |
|                                                                                     |  | Usable Length<br><b>95, 105, 115cm</b>    | Distal Hydrophilic Coating               |                     |
|                                                                                     |                                                                                    |                                           | BER<br><b>11.5cm</b>                     | SIM2<br><b>30cm</b> |
|                                                                                     | French Designation (FD)<br><b>7F</b>                                               | Outer Diameter (OD)<br><b>0.096in</b>     | Internal Diameter (ID)<br><b>0.087in</b> |                     |
|  | Usable Length<br><b>95, 105cm</b>                                                  | Distal Hydrophilic Coating<br><b>25cm</b> |                                          |                     |

COMPONENTS

Peel-Away  
Introducer

Custom BERN  
or SIM2

Dilator



UNPARALLELED  
TRACKING AND SUPPORT



PATENTED INNOVATIVE  
TECHNOLOGY



SWITCH BETWEEN MODES  
ON-DEMAND

Our patented technology places full **unparalleled control** in the palm of your hand. The **first-of-its-kind design** allows operators to seamlessly adapt the catheter to strike the **ideal balance between trackability and support** on demand as it navigates through complex neurovascular anatomy.

More trackable

More control

More deliverable

Better navigability

Stronger support

# Powerful. Reliable. Efficient.

We create bold solutions to challenge AFib. You rely on technology bold enough not to compromise between efficiency and accuracy. With this technology that boldly redefines the EP Lab, you can elevate your practice and meet today's challenges.

Designed to revolutionize high-density cardiac mapping, our NEW investigational EP System and Technology.

Our product is a next generation device currently undergoing clinical trials, is currently being studied in clinical studies for safety and effectiveness.

Product has an approved IDE, investigational device exemption

Our investigational device gained IDE approval from the FDA and is currently enrolling patients.

CAUTION: Investigational device. Limited by Federal (or United States) law to investigational use.

# Powerful. Reliable. Efficient.

We create bold solutions to challenge AFib. You rely on technology bold enough not to compromise between efficiency and accuracy. With this technology that boldly redefines the EP Lab, you can elevate your practice and meet today's challenges.

Designed to revolutionize high-density cardiac mapping, our NEW investigational EP System and Technology.

Our product is a next generation device currently pending 510(k) clearance\*.

Product has an approved IDE, investigational device exemption

Our investigational device gained IDE approval from the FDA and is currently enrolling patients.

# Clinical case cines / radiological images

